



BITE REGISTRATION

Mousse-like vinyl polysiloxane derivative material in a high performance cartridge dispensing system. Hard setting material with excellent tear strength.

FEATURES

Nonflowing, thixotropic

Available in three different set times

Hard set

BENEFITS

No drip, sag or slump, stays in place on occlusal tooth surfaces

Reduces patient discomfort and speeds use of chair time

Prevents distortion in plastic, disposable bite trays

	WORK TIME*	INTRAORAL SET TIME	SHORE "A"
Regular Set	2:00	1:00	—
Fast Set	0:30	1:00	—
Super Fast Set	0:30	0:30	—
Rigid Bite			
Fast Set	0:30	1:00	90+

*Includes mixing time.

Contains: 2 x 50 ml cartridges and 6 x 6.5 mm mixing tips.

CODE	DESCRIPTION
Unflavored	
Q-9001HP	Regular Set
Q-9002HP	Fast Set
Q-9000HP	Super Fast Set
Flavored - Vanilla Crème	
Q-9003HP	Fast Set
Q-9004HP	Super Fast Set
Rigid Bite	
Q-9005HP	Rigid Bite - Fast Set, Unflavored



Distributed by:



Quala
DENTAL PRODUCTS
www.quala.com

#QSS035R2