



# CLEAR BITE REGISTRATION HIGH PERFORMANCE

Clear vinyl polysiloxane derivative material in a high performance cartridge dispensing system. Clear Bite Registration can be used as a matrix material for temporaries.

## FEATURES

Clear “see-through” material

Non-flowing, thixotropic

## BENEFITS

Allows light to pass through in curing of temporaries, makes any bubbles or voids visible

No drip, sag or slump, stays in place on occlusal tooth surfaces

	WORK TIME*	INTRAORAL SET TIME	SHORE “A”
Fast Set	1:00	1:00	60

\* Includes mixing time.

*Contains: 2 x 50 ml cartridges and 6 x 6.5 mm mixing tips.*

CODE	DESCRIPTION
------	-------------

**Unflavored**

<b>Q-9006HP</b>	Fast Set
-----------------	----------



Distributed by:



**Quala**  
DENTAL PRODUCTS  
www.quala.com

#QSS077R2