



EXTRA SOFT EARLOOP FACE MASK

Extra soft material for sensitive skin wearers. Cellulose inner layer is made of hypoallergenic material that is void of dyes, inks and chemicals. The inner layer, which does not lint, tear or shred is combined with the spunbond fluid-resistant outer layer to provide the ultimate protection. Mask meets ASTM Level 3 specifications for maximum fluid-resistant protection with a PFE rating of greater than 99% at 0.1 microns. Latex-Free.



CODE	DESCRIPTION	COLOR	QTY
GCUWHA	Extra Soft Earloop Mask	White	50/box

Distributed by:

