

Quala Dental Products

Instruments



Made in the USA

Titanium Implant Scalers & Curettes

Used on implant abutments to remove calculus cleanly and easily. Titanium tips are thin to access tight interproximal spaces and easy to adapt along the gingival margin. Stainless steel construction provides maximum durability.

FEATURES

3/8" diameter light weight handle

Solid, Class 5 Titanium tips

Thin blades

Blades never require sharpening

BENEFITS

Relaxed, ergonomic grip

Won't scratch implant abutments, glides easily over surface

Provides better access to tight spaces than plastic implant scalers

Implant abutments won't dull instrument

CODE	DESCRIPTION
SIH6H7/6Q	Scaler, Anterior Offset Sickle 6/7, #6 Handle
SI204S/6Q	Scaler, Posterior Sickle 204S, #6 Handle
SIBH5-6/6Q	Curette, Barnhart 5-6, #6 Handle

Distributed by: