



INTRAORAL DENTAL FILM

A high contrast film, providing greater resolution and detail, available in both D Speed and F Speed for reduced exposure. The flexible vinyl film packets have easy-opening tabs, with ultra soft edges for increased patient comfort. Films can be processed manually or automatically. 24-month shelf life.

FEATURES

High contrast film

Soft, supple vinyl packets

Easy opening tabs

BENEFITS

Greater resolution and detail

Patient comfort

For effortless film removal

CODE	DESCRIPTION	QTY
18200-QU	QD 58 Adult D Speed	150/box
18310-QU	QF 58 Adult F Speed	150/box
18201-QU	QD 57 Adult D Speed, Double Packs	130/box
18410-QU	QF 57 Adult F Speed, Double Packs	130/box
18235-QU	QD 54 Child D Speed	100/box
18635-QU	QF 54 Child F Speed	100/box



WARNING: Intraoral Dental Film – This product contains lead, which is known to the state of CA to cause cancer, birth defects, or other reproductive harm. For more information go to www.P65Warnings.ca.gov.

Distributed by:



Quala
DENTAL PRODUCTS
www.quala.com

#QSS073R3