



# Carbide Burs

## Operative Carbides

	Round					Inverted Cone				Straight Flat End Fissure	
<b>Shape No.</b>	1/2	2	4	6	8	33 1/2	34	35	37	56	57
<b>Size (1/10 mm)</b>	006	010	014	018	023	006	008	010	014	009	010
<b>Length (mm)</b>	-	-	-	-	-	0.4	0.6	0.7	0.8	3.2	3.7
<b>Friction Grip</b>											
10-Pack	Q2060001	Q2060050	Q2060086	Q2060104	Q2060154	Q2060072	Q2060079	Q2060080	Q2060082	Q2060096	Q2060099
50-Pack	-	Q2050050	Q2050086	Q2050154	Q2050154	-	-	Q2050080	-	Q2050096	Q2050099
<b>Latch</b>											
10-Pack	-	-	-	-	-	-	-	-	-	-	-
50-Pack	-	Q2060339	Q2060353	Q2060363	Q2060375	-	-	-	-	-	-

## Operative Carbides

	Taper Flat End Fissure	Amalgam Prep	Pear				Straight Flat End X-Cut Fissure			Taper Flat End X-Cut Fissure
<b>Shape No.</b>	170	245	329	330	331	332	556	557	558	701
<b>Size (1/10 mm)</b>	010	009	006	008	010	012	009	010	012	012
<b>Length (mm)</b>	3.7	2.7	1.6	1.6	1.6	1.6	3.2	3.7	3.7	3.7
<b>Friction Grip</b>										
10-Pack	Q2060026	Q2060062	Q2060071	Q2060073	Q2060074	Q2060076	Q2060089	Q2060091	Q2060093	Q2060117
50-Pack	-	Q2050062	-	Q2050073	Q2050074	-	Q2050089	Q2050091	Q2050093	-

## Operative Carbides

	Straight Round End Fissure		Straight Round End X-Cut Fissure		Trimming & Finishing (Gold-Plated)					
					Egg			Needle		
<b>Shape No.</b>	1157	1158	1557	1558	7404	7406	7408	9406	7901	7902
<b>Size (1/10 mm)</b>	010	012	010	012	014	018	023	018	009	010
<b>Length (mm)</b>	3.7	3.7	3.7	3.7	2.8 (12 blade)	3.4 (12 blade)	3.6 (12 blade)	3.4 (30 blade)	2.2 (12 blade)	3.4 (12 blade)
<b>Friction Grip</b>										
10-Pack	Q2060007	Q2060008	Q2060018	Q2060020	Q2060134	Q2060135	Q2060136	Q2060161	Q2060151	Q2060152
50-Pack	Q2050007	Q2050008	-	-	-	-	-	-	-	-



**Qualala**<sup>®</sup>  
**DENTAL PRODUCTS**