



ALGINATE SUBSTITUTE

For use with preliminary orthodontic models, provisionals, impressions for temporaries and many other applications. Quala Alginate Substitute material produces accurate impressions every time. Impressions can be placed in disinfectant solution without damage.

FEATURES

- Pre-mixed material
- Excellent dimensional stability
- Pleasant mint flavor

BENEFITS

- No messy mixing of powder and liquid
- Multiple pours can be made without loss of accuracy
- Increases patient acceptance, helps to reduce gag reflexes during impressions

	WORK TIME	INTRAORAL SET TIME	SHORE "A"
Fast Set	1:00 minute	1:00 minute	50

CODE DESCRIPTION

Alginate Substitute

- Q589510** 50 mL Fast Set, Mint flavor
Contains: 8 x 50 mL cartridges
- Q589512** 380 mL Fast Set, Mint flavor
Contains: 2 x 380 mL cartridges

CODE DESCRIPTION QTY

Alginate Substitute

- Q-8020** Dynamic Mixing Tips, Yellow 50/bag
- Q-8021** Bayonet Ring, Yellow 2/bag



Distributed by:

