



ASPIRATOR TIPS

Rigid plastic construction provides greater cheek retraction. Pre-bent 30° angle. Available in both a 1/4" and 1/8" diameter for macro-suctioning or micro-aspiration.

FEATURES


Manufactured of rigid polyethylene

Pre-bent 30° angle

BENEFITS

Provides greater cheek retraction. Dispose after single use, or alternately, tips are autoclavable up to 275° F.

Greater visual access



| CODE | DESCRIPTION | COLOR | QTY |
|----------|---------------|-------|--------|
| 403201-6 | 1/8" Diameter | White | 25/bag |
| 403202-6 | 1/4" Diameter | Green | 25/bag |

Distributed by:

