



BIOLOGICAL MONITORING SYSTEM

Convenient, reliable system for monitoring effectiveness of sterilization. Bacterial spore test strips may be used to test all sterilization processes: Steam, Chemical Vapor (Chemiclave), Dry Heat and Ethylene Oxide.

FEATURES

Mail-in service

Third party biological indicator effectiveness documentation

Complimentary reminder program

Failure results communicated immediately

BENEFITS

Confidential test results within 24 hours after receipt at laboratory with ability to check results online

Provides proof of verification results

Encourages routine testing compliance, as suggested by CDC

Office can take immediate action

CODE	DESCRIPTION	SIZE
Monthly		
Q602412	Contains: Two test strips and one control strip; Pre-paid postage included	12 tests/box
Weekly		
Q602448	Contains: Two test strips and one control strip; Pre-paid postage included	48 tests/box
Value		
Q602452	Contains: One test strip and one control strip; Postage not included	52 tests/box

Made in the USA.

Distributed by:

