



CARBIDE BURS

Carbide Burs provide superior cutting performance of restorative materials. The shank of each bur is nickel or gold plated to protect against corrosion and to provide a firm grip in the handpiece. Variety of shapes and sizes allows for more flexibility in performing procedures. Available in 10-packs and 50-packs based on select styles. Autoclavable.

Operative Carbides	Round					Inverted Cone				Straight Flat End Fissure	
Shape No.	1/2	2	4	6	8	33 1/2	34	35	37	56	57
Size (1/10 mm)	006	010	014	018	023	006	008	010	014	009	010
Length (mm)	-	-	-	-	-	0.4	0.6	0.7	0.8	3.2	3.7
Friction Grip											
10-Pack	Q2060001	Q2060050	Q2060086	Q2060104	Q2060154	Q2060072	Q2060079	Q2060080	Q2060082	Q2060096	Q2060099
50-Pack	-	Q2050050	Q2050086	Q2050154	Q2050154	-	-	Q2050080	-	Q2050096	Q2050099
Latch											
10-Pack	-	Q2060339	Q2060353	Q2060363	Q2060375	-	-	-	-	-	-

Operative Carbides	Taper Flat End Fissure	Amalgam Prep	Pear				Straight Flat End X-Cut Fissure			Taper Flat End X-Cut Fissure
Shape No.	170	245	329	330	331	332	556	557	558	701
Size (1/10 mm)	010	009	006	008	010	012	009	010	012	012
Length (mm)	3.7	2.7	1.6	1.6	1.6	1.6	3.2	3.7	3.7	3.7
Friction Grip										
10-Pack	Q2060026	Q2060062	Q2060071	Q2060073	Q2060074	Q2060076	Q2060089	Q2060091	Q2060093	Q2060117
50-Pack	-	Q2050062	-	Q2050073	Q2050074	-	Q2050089	Q2050091	Q2050093	-

Operative Carbides	Straight Round End Fissure		Straight Round End X-Cut Fissure		Trimming & Finishing (Gold-Plated)				Needle	
Shape No.	1157	1158	1557	1558	7404	7406	7408	9406	7901	7902
Size (1/10 mm)	010	012	010	012	014	018	023	018	009	010
Length (mm)	3.7	3.7	3.7	3.7	2.8 (12 blade)	3.4 (12 blade)	3.6 (12 blade)	3.4 (30 blade)	2.2 (12 blade)	3.4 (12 blade)
Friction Grip										
10-Pack	Q2060007	Q2060008	Q2060018	Q2060020	Q2060134	Q2060135	Q2060136	Q2060161	Q2060151	Q2060152
50-Pack	Q2050007	Q2050008	-	-	-	-	-	-	-	-

WARNING: This product contains nickel, a chemical known to the State of California to cause cancer. For more information go to www.P65warnings.ca.gov

Distributed by:

