



DENTAL AMALGAM ALLOY

Dispersed Phase High Copper

Dispersed Phase High Copper Dental Amalgam Alloy creates a strong, long-lasting, alloy ideal for posterior restoration. The unique high copper formula offers superior handling characteristics, early high-compressive strength, with a smooth carve. Approximate metal composition: 41.5% silver, 30.5% tin, 28% copper.

FEATURES

- High early strength
- Self-activating capsule
- Color-coded capsules

BENEFITS

- Reduces early breakage
- Requires no pinching or twisting to open. One-step trituration without the use of an activator tool
- Easy differentiation between spills

	WORK TIME	SET TIME
Fast Set	2:00 minute	5-6 minutes

CODE	DESCRIPTION
Fast Set, 50/pkg	
Q055022	1 Spill, Green capsules, 400 mg
Q055026	2 Spill, Blue capsules, 600 mg
Q055030	3 Spill, Yellow capsules, 800 mg
Fast Set, Bulk 500/pkg	
Q055003	1 Spill, Green capsules, 400 mg
Q055007	2 Spill, Blue capsules, 600 mg

WARNING: This product contains mercury. Can cause birth defects and other reproductive harm. www.P65Warnings.ca.gov.

Distributed by:

