



DISPOSABLE DENTAL INJECTION NEEDLES

Pre-threaded plastic hub needles. Bevel indicator on the hub for quick orientation of the needle. Color-coded caps indicate gauge size.

FEATURES

Pre-threaded Plastic Hubs

Bevel Indicator

Siliconized Cannula

Color-coded Caps

BENEFITS

Easy assembly with less risk of damage to syringes

For quick orientation of the needle for different injections

Ensures a smooth insertion

Easy identification of gauge size



CODE	DESCRIPTION	COLOR	QTY
DP47825	25 Gauge Short 1"	Pink	100/box
DP40825	25 Gauge Long 1 3/8"	Red	100/box
DP47827	27 Gauge Short 1"	Orange	100/box
DP47827	27 Gauge Long 1 3/8"	Yellow	100/box
DP47830	30 Gauge Short 1"	Blue	100/box
DP42830	30 Gauge X-Short 3/8"	Purple	100/box

Distributed by:

