



## EARLOOP MASKS

Procedural mask with soft white facial tissue inner layer, fluid resistant outerfacing, and outside earloop attachment to eliminate potential irritation. Particle filtration efficiency of 99.5% at 1 micron. Latex-free.

### FEATURES

Filtration

Fluid resistant outerfacing

Outside earloop attachment

### BENEFITS

99.5% efficiency at 1 micron

Resists fluid penetration to wearer

Eliminates irritation

CODE	DESCRIPTION	QTY
GCI PKA	Pink Earloop	50/box
GCI BLA	Blue Earloop	50/box
GCI LVA	Lavender Earloop	50/box



Distributed by:

