



## ENZYME ULTRASONIC TABLETS

For general purpose ultrasonic cleaning, evacuation system cleaning and instrument presoak. Effective, fast removal of blood, food particles, compounds and other soils. Two tablets make one gallon of solution.  
Clean mint scent.

### FEATURES

Compact and convenient,  
2 tablets mix with water  
to make 1 gallon of solution

Drop tablets directly into  
ultrasonic cleaner

Clean mint scent

Contains fast-acting  
protease enzymes

Nonionic, low foaming,  
biodegradable

Also effective as  
Evacuation Cleaner

### BENEFITS

Saves counter and storage space;  
No lifting of heavy liquids

No messy mixing of liquid  
concentrates; easy to use

Pleasant to work with;  
minimizes any odor from tank

Effective, fast removal of blood,  
particles, compounds and  
other soils which are normally  
resistant to detergents alone

Easy on instruments and  
the environment

One product for multiple needs



CODE	QTY
50036380	64 tablets/box

Distributed by:

