



GENERAL PURPOSE ULTRASONIC SOLUTION

Low foaming solution effectively removes pumice, rouge, buffing compounds, food particles, foreign matters and soils of all types from instruments and appliances in all ultrasonic cleaners. Dilute 1 part solution to 10 parts water. Use directly in ultrasonic tank. Solution color: Blue.

FEATURES

General purpose formula with isopropanol

Nonionic, nonsudsing

10:1 concentrate

BENEFITS

Effective soil remover in all ultrasonic cleaners

Low foaming action, rapid wetting; eliminates messy cleaning

Greater economy; easier storage



CODE	DESCRIPTION
50036360	General Purpose Ultrasonic Solution , 1 gallon

Distributed by:

