



PLASTER & STONE REMOVER

For removal of plaster and stone from bridges, crowns, dentures, inlays, metal impression trays, spatulas, etc. Use at full strength in glass beaker placed inside ultrasonic cleaner, 5-10 minutes only. Solution color: Pink.

FEATURES

Efficient, effective, safe

Contains rust inhibitor

Nonionic

BENEFITS

Eliminates hand scrubbing of instruments; saves time; protects staff from potential contamination from instruments

Protects instruments from corrosion, extends instrument life

Low foam level for better, faster, more complete cleaning



CODE	DESCRIPTION
50036365	Plaster & Stone Remover, 1 gallon

Distributed by:

