



SELF-SEALING STERILIZATION POUCHES

Dual Indicator sterilization pouches feature process indicators that are printed internally and externally on each pouch and change color during sterilization to verify pouches have been processed. These process indicators only measure temperature. Pouches are made from medical grade paper and feature blue tint film, solid seal and thumb notch. Latex-free.

FEATURES

Wide solid seal

Dual process indicators for steam or EtO

See-through four-ply film

Tack seal over perforation

Edge overlap

Steam absorbent paper

BENEFITS

Protects against separation, puncture and instrument contamination

Pink arrows will turn brown when processed by steam. Blue arrows will turn gold/brown when EtO processed.

Lets you identify instrument set-up; helps reveal punctures

Prevents separation of sealing end, and prevents pouch from being compromised

Allows easy pouch opening to access contents

Allows steam through, yet remains strong even when wet

CODE	DESCRIPTION	CODE	DESCRIPTION
200/box		Bulk 500/box	
AD2440	2 ¼" x 4"	LAD2790	2 ¾" x 9"
AD2790	2 ¾" x 9"	LAD3590	3 ½" x 9"
AD3554	3 ½" x 5 ¼"	LAD5410	5 ¼" x 10"
AD3590	3 ½" x 9"		
AD5410	5 ¼" x 10"		



Distributed by:

