



STERILIZATION WRAP

A soft, flexible material to wrap instruments for autoclave sterilization.
Blue, autoclavable, disposable.

FEATURES

Produced from natural cellulose fibers

Compatible with various methods of sterilization

BENEFITS

Provides the highest barrier to airborne and waterborne bacterial contamination

Use with steam, EtO, gamma irradiation, and E-beam



CODE	DESCRIPTION	QTY
Q602330	15" x 15"	100 sheets/bag
Q602335	20" x 20"	100 sheets/bag
Q602340	24" x 24"	100 sheets/bag

Distributed by: