



## TARTAR & STAIN REMOVER

Removes calculus, tartar, tobacco, and food stains from dentures, partials, crowns and bridges. Use at full strength in glass beaker placed inside ultrasonic cleaner, 5-10 minutes only. Solution color: Clear.

### FEATURES

Removes calculus, tartar, food stains and cement from dentures, partials, crowns, bridges

Use at full strength

### BENEFITS

Faster cleanup time; less handling

No dilution, no mixing necessary



CODE	DESCRIPTION
50036370	Tartar & Stain Remover, 1 gallon

Distributed by: