



## TEMPORARY CEMENT REMOVER

Formulated for removal of temporary cement from bridges, crowns, inlays, dentures, orthodontic appliances, instruments, and other items. Use at full strength in glass beaker placed inside ultrasonic cleaner, 5-10 minutes only. Solution color: Yellow.

### FEATURES

Noncorrosive

For removal of temporary cements including eugenol based

### BENEFITS

Protects instruments, extends instrument life

Minimizes cleanup time



CODE	DESCRIPTION
50036375	Temporary Cement Remover, 1 gallon

Distributed by:

