



TRIPLE ENZYMATIC EVACUATION CLEANER

Non-foaming formula cleans evacuation systems and deodorizes as it cleans. Triple enzyme formula (Protease; Amylase; Lipase) dissolves organic fluids and debris to restore maximum performance to system and hoses. Neutral pH formula leaves a clean, fresh scent.

FEATURES

Economical concentrate

Environmentally friendly, neutral pH

Easy dispensing

Safe to use every day

BENEFITS

One ounce mixes with warm water to make 32 oz. of solution

Safe for use with amalgam separators

Solution squeezes easily into pre-measured reservoir

Daily use improves cleaning efficiency



CODE	DESCRIPTION
Q182605	1 Quart (32 oz.) Concentrate

Distributed by:



Quala
DENTAL PRODUCTS
www.quala.com

#QSS081R2