



ULTRASONIC CLEANER WITH TIMER

Quala Ultrasonic Cleaners feature self-contained generators in a variety of sizes and capacities to provide efficiency and dependability for all your cleaning applications. Constructed with a durable stainless steel tank, with vinyl-clad steel and aluminum housing. Equipped with a 60-minute timer. 2½ year warranty on all cleaners.

FEATURES

Compact size

Swivel drain

Cover included

BENEFITS

Perfect for countertop use

Multi-positional outlet for easy removal of solution

Reduces the possibility of water splashing to work area



CODE	DESCRIPTION
5002 Cleaner with Timer	
410	Internal tank dimensions: 9 $\frac{3}{8}$ "x 5 $\frac{3}{8}$ "x 4" (23.8 x 13.7 x 10.2 cm) Capacity: 3 $\frac{3}{4}$ quarts (3.2 liters)
5002 Accessories	
423	Accessory Kit: Includes one solid side basket, positioning cover, two 600 mL beakers complete and one bur tray.
1351	Stainless Steel Basket for 5002

Distributed by: