



## SUTURES

Precision drilled end needle ensures smooth transition to suture material and maximizes secure attachment of suture to needle. Needles are sharpened to a smooth, precise cutting edge to penetrate cleanly and stay sharp for repeated passage through tissue. Double sterile pouch allows sterile transfer onto a sterile field. 12 per box.

### FEATURES

Unique needle design and manufacture

See-through packaging



Industry standard label colors

### BENEFITS

Needles stay sharp for repeated passage through tissue

Saves money by reducing selection errors and wasted materials

Simplifies selection, enhances efficiency





Item#	Needle Code	Needle Description	Suture Size	Suture Length
-------	-------------	--------------------	-------------	---------------

### Chromic Gut, Reverse Cutting

P00687	C-3	3/8 circle – 13 mm	5/0	18" (45 cm)
P00558	C-6	3/8 circle – 19 mm	4/0	18" (45 cm)
P00525	C-6	3/8 circle – 19 mm	3/0	18" (45 cm)
P00635	C-6	3/8 circle – 19 mm	4/0	27" (67.5 cm)
P00636	C-6	3/8 circle – 19 mm	3/0	27" (67.5 cm)

### Nylon Black Monofilament, Reverse Cutting

P00662	C-6	3/8 circle – 19 mm	4/0	18" (45 cm)
--------	-----	--------------------	-----	-------------

### Plain Gut, Reverse Cutting

P00592	C-6	3/8 circle – 19 mm	3/0	18" (45 cm)
P00071	C-6	3/8 circle – 19 mm	4/0	18" (45 cm)
P00822	C-6	3/8 circle – 19 mm	3/0	27" (67.5 cm)
P00810	C-7	3/8 circle – 24 mm	3/0	27" (67.5 cm)
P00065	C-31	1/2 circle – 24 mm	3/0	18" (45 cm)

### Silk Black Braided, Reverse Cutting

P00066	C-6	3/8 circle – 19 mm	4/0	18" (45 cm)
P00072	C-6	3/8 circle – 19 mm	3/0	18" (45 cm)
P00062	C-26	1/2 circle – 15 mm	3/0	18" (45 cm)
P00064	C-31	1/2 circle – 24 mm	3/0	18" (45 cm)

Made in Mexico.

Distributed by:

