



# VINYL POLYSILOXANE (VPS)

Super hydrophilic formulation offers excellent dimensional stability and tear strength for crisply detailed impressions every time.

## FEATURES

State-of-the-art hydrophilic technology

Uniform consistency

High dimensional stability

Pleasantly neutral taste and odor

## BENEFITS

Provides great “wetability”; it won’t “run away” from prepared impression area

Assures crisply detailed impressions

Maintains the accuracy of the impression over time, even after disinfection by immersion

Increases patient acceptance, helping to reduce gag reflexes during impressions



|             | WORK TIME* | INTRAORAL SET TIME |
|-------------|------------|--------------------|
| Regular Set | 2:30       | 1:00               |
| Fast Set    | 1:00       | 1:00               |

\* Includes mixing time.

Contains: 4 x 50 ml cartridges. Tips sold separately.

| CODE | DESCRIPTION |
|------|-------------|
|------|-------------|

### Light Body (Yellow)

**32823C** Regular Set, Unflavored

**32822C** Fast Set, Unflavored

**32821C** Fast Set, Vanilla Crème

### Medium Body/Monophase (Purple)

**32826C** Regular Set, Unflavored

**32825C** Fast Set, Unflavored

**32824C** Fast Set, Vanilla Crème

### Heavy Body (Rose)

**32829C** Regular Set, Unflavored

**32828C** Fast Set, Unflavored

**32827C** Fast Set, Vanilla Crème



Distributed by:

