



# ENZYMATIC DETERGENT CONCENTRATE

Enzymatic detergent, presoak and evacuation system cleaner removes blood, protein, tissue and soil from instruments, glassware, equipment, linens, clothing and evacuators. Pleasant spearmint scent.

### FEATURES

Concentrated

Advanced corrosion protection

Rapid cleaning

Bioburden reduction

Multipurpose

Environmentally friendly

### BENEFITS

Mix 1-2 oz. of solution to 1 gallon of water

Unique anticorrosion formulation will not damage instruments or devices

Enzymatic detergent quickly and easily removes blood and soil in minutes without scrubbing

Unique dual-enzyme formula assists in reducing Gram-negative, Gram-positive bacteria and viruses

This product can be used as a presoak, manual detergent, ultrasonic detergent, automatic washer detergent and linen presoak. Use twice weekly to keep evacuators clean and odor free

Formula is biodegradable and phosphate free



CODE

SIZE

18-555

1 gallon

Distributed by:



**Quala**  
DENTAL PRODUCTS  
[www.quala.com](http://www.quala.com)

#QSS052R2