Quala®
DENTAL PRODUCTS

DISPOSABLES

NON-WOVEN SPONGES

Rayon/polyester blend for a softer texture with increased volume and absorbency.

FEATURES

70% rayon, 30% polyester

BENEFITS

Softer texture, increases the volume and absorbency



CODE	DESCRIPTION	QTY
ENCNWA	2” x 2”, 4-Ply (50g)	200 /sleeve 25 sleeves/case (5000)
ENC4NWA	4” x 4”, 4-Ply (40g)	200 /sleeve 10 sleeves/case (2000)

Distributed by:

Quala®
DENTAL PRODUCTS
www.quala.com
#QSS024R3

© 2025 NDC, Inc.