Quala®
DENTAL PRODUCTS

DISPOSABLES

COTTON-FILLED SPONGES

Manufactured with cotton gauze made from pure virgin white cotton with a 100% cotton fiber filling for maximum absorbency. Non-sterile.

FEATURES	BENEFITS
100% pure virgin white cotton	No impurities
100% cotton fiber filling	Maximum absorbency
"Finer" gauze weave	Enhances rate and volume of absorbency
Superior material and construction	Minimizes "dust"
USP Type VII gauze	Tighter 20 x 12 mesh



CODE	DESCRIPTION	QTY
ENCCA	2" x 2" Cotton-Filled Sponges	200 /sleeve 25 sleeves/case (5000)

Distributed by:

Quala®
DENTAL PRODUCTS
www.quala.com
#QSS023R2

© 2018 NDC, Inc.