DIAMONDS

Variety of shapes and grits for more flexibility in performing procedures. Diamond particles are uniform in coverage for optimal cutting performance, less heat and a longer lasting diamond. A concentric stainless steel shank means minimal runout, greater stability at higher speeds. Available in Single-Use with grits ranging from Fine to Super Coarse..

Single-Use Diamonds

Available in Packs of 10 & 25

	Occlusal Reduction/Finishing						Flat End Cylinder						
	Football		Large Head Football		Large Head								
Shape No.	368	368	368	379	811	811	835	835	835	835	836	837	S835
Size (1/10 mm)	016	018	023	023	033	037	008	009	010	012	012	012	010
Length (mm)	3.0	5.0	5.0	4.5	4.2	7.0	4.0	3.0	4.0	4.0	5.0	7.0	4.0
Extra Fine			0032										
Fine	0026	0029	0002	0033									
Medium									0040				0045
Coarse	0027	0030	0001	0004	0034	0036	0038	0039	0006	0007	0041	0043	0046
Super Coarse	0028	0031	0003	0005	0035	0037					0042	0044	

Single-Use Diamonds

Available in Packs of 10 & 25

	Flat End Taper							Round End Taper											
							Large Head						Short Shank						Large Head
Shape No.	847	847	847	847	848	848	848	850	850	850	855	855	855	855	856	856	856	856	856
Size (1/10 mm)	012	014	016	018	016	018	023	012	014	018	012	016	016	018	014	016	018	021	023
Length (mm)	9.0	9.0	9.0	8.0	11.0	10.0	10.0	11.0	10.0	10.0	7.0	7.0	7.0	7.0	8.0	8.0	9.0	8.0	7.0
Extra Fine								0062											
Fine	0047							0063	0065	0067				0075		0079	0081		
Medium		0048	0050	0053		0056	0059	0064	0066	0068		0071	0072	0076	0078		0082		
Coarse	0008	0049	0051	0054	0055	0057	0060	0010	0011	0069	0070	0073	0014	0077	0017	0012	0024	0083	0085
Super Coarse			0052	0009		0058	0061			0013		0074	0015	0016		0080	0025	0084	0086

Single-Use Diamonds

Available in Packs of 10 & 25

	Needle						Flame									Beveled
	Large Head					Short Shank										
Shape No.	852	858	858	859	889	858	862	862	863	863	863	863	879	878K	879K	878
Size (1/10 mm)	023	008	014	014	008	008	012	016	012	014	016	018	018	012	012	012
Length (mm)	10.0	3.0	8.0	10.0	3.0	3.0	8.0	8.0	10.0	8.0	10.0	9.0	9.0	8.0	10.0	8.0
Extra Fine		0089			0096	0090	0099	0103	0107							
Fine	0087		0091	0018	0097		0100	0104	0108	0110	0112	0114	0116			0123
Medium			0092	0094	0098		0101	0105	0109					0119	0121	0124
Coarse	0088		0093	0095			0102	0106	0019	0020	0113	0115	0117	0120	0122	0021
Super Coarse										0111			0118			0022

Diamonds Ordering Information

When ordering diamonds, use the following pkg quantity prefix from the list at the right in addition to the item number from the charts above for the total item code. For example, to order a 10-pack of 835-008 flat end cylinder coarse single-use diamonds, use part number Q2080038 (Q208 + 0038).

Single-Use Diamonds
Q208 = 10/pkg
Q209 = 25/pkg

Quala® DENTAL PRODUCTS

DIAMONDS

Single-Use Diamonds	Cylinder		Bulk Reducer	Curettage				Double Cone	
Available in Packs of 10 & 25									
Shape No.	**878**	**879**	**6055**	**878K**	**878K**	**878K**	**879K**	**806**	**813**
Size (1/10 mm)	014	012	018	014	016	018	018	012	014
Length (mm)	8.0	10.0	8.0	9.0	9.0	8.0	11.0	2.5	2.0
Extra Fine									
Fine									
Medium	0125	0127		0130	0132	0134	0136	0138	0139
Coarse	0126	0128	0129	0131	0133	0135	0137		
Super Coarse			0023						

Diamonds Ordering Information

When ordering diamonds, use the following pkg quantity prefix from the list at the right in addition to the item number from the chart below for the total item code. For example, to order a 10 pack of 835-008 flat end cylinder coarse single-use diamonds, use part number Q2080038 (Q208 + 0038).

Single-Use Diamonds
Q208 = 10/pkg
Q209 = 25/pkg

Turbo Diamonds
Q208 = 10/pkg

Single-Use Diamonds	Inverted Cone			Modified Shoulder		Mosquito	Pear		Round			Round End Cylinder				Wheel		Depth Cutter
											Large Head					Large Head		
Available in Packs of 10 & 25																		
Shape No.	**805**	**805**	**805**	**837KR**	**850KR**	**392**	**830**	**830**	**801**	**801**	**801**	**838**	**880**	**880**	**881**	**909**	**909**	**834**
Size (1/10 mm)	010	014	016	012	018	016	010	014	010	014	023	012	010	012	014	035	043	021
Length (mm)	1.0	2.0	2.0	8.0	9.0	5.0	2.2	3.0	1.0	1.4	2.3	4.0	7.0	7.0	8.0	1.1	1.3	6.0
Extra Fine										0156								
Fine						0147												
Medium	0140	0142				0148	0149	0151	0154	0157	0166							0173
Coarse	0141	0143	0144	0145	0146		0150	0152	0155	0158		0159	0160	0162	0164	0169	0171	
Super Coarse								0153					0161	0163		0170	0172	

Turbo Diamonds	Occlusal Reduction/Finishing	Flat End Taper		Round End Taper					Flame	Modified Shoulder		Round End Cylinder
Available in Packs of 10												
Shape No.	**368T**	**847T**	**847T**	**850T**	**850T**	**855T**	**855T**	**856T**	**863T**	**837KR**	**850KRT**	**880T**
Size (1/10 mm)	023	016	018	014	018	016	018	018	014	012	018	010
Length (mm)	5.0	9.0	8.0	10.0	10.0	7.0	7.0	9.0	8.0	8.0	9.0	7.0
Item #	Q2080174	Q2080175	Q2080176	Q2080177	Q2080178	Q2080179	Q2080180	Q2080181	Q2080182	Q2080183	Q2080184	Q2080185

Distributed by:

Quala®
DENTAL PRODUCTS
www.quala.com
#QSS070R3