DIAMONDS

Single-Use Diamonds

Available in Packs of 10 & 25

Occlusal Reduction/Finishing



	Football		Large Head Football		Large Head	
Shape No.	368	368	368	379	811	811
Size (1/10 mm)	016	018	023	023	033	037
Length (mm)	3.0	5.0	5.0	4.5	4.2	7.0
Extra Fine			0032			
Fine	0026	0029	0002	0033		
Medium						
Coarse	0027	0030	0001	0004	0034	0036
Super Coarse	0028	0031	0003	0005	0035	0037

Flat End Cylinder



Shape No.	835	835	835	835	836	837	S835
Size (1/10 mm)	008	009	010	012	012	012	010
Length (mm)	4.0	3.0	4.0	4.0	5.0	7.0	4.0
Extra Fine							
Fine							
Medium			0040				0045
Coarse	0038	0039	0006	0007	0041	0043	0046
Super Coarse					0042	0044	

Flat End Taper



							Large Head
Shape No.	847	847	847	847	848	848	848
Size (1/10 mm)	012	014	016	018	016	018	023
Length (mm)	9.0	9.0	9.0	8.0	11.0	10.0	10.0
Extra Fine							
Fine	0047						
Medium		0048	0050	0053		0056	0059
Coarse	0008	0049	0051	0054	0055	0057	0060
Super Coarse			0052	0009		0058	0061

Round End Taper



						Short Shank						Large Head
Shape No.	850	850	850	855	855	855	855	856	856	856	856	856
Size (1/10 mm)	012	014	018	012	016	016	018	014	016	018	021	023
Length (mm)	11.0	10.0	10.0	7.0	7.0	7.0	7.0	8.0	8.0	9.0	8.0	7.0
Extra Fine	0062											
Fine	0063	0065	0067				0075		0079	0081		
Medium	0064	0066	0068		0071	0072	0076	0078		0082		
Coarse	0010	0011	0069	0070	0073	0014	0077	0017	0012	0024	0083	0085
Super Coarse			0013		0074	0015	0016		0080	0025	0084	0086

Single-Use Diamonds

Available in Packs of 10 & 25

Needle



	Large Head					Short Shank
Shape No.	852	858	858	859	889	858
Size (1/10 mm)	023	008	014	014	008	008
Length (mm)	10.0	3.0	8.0	10.0	3.0	3.0
Extra Fine		0089			0096	0090
Fine	0087		0091	0018	0097	
Medium			0092	0094	0098	
Coarse	0088		0093	0095		
Super Coarse						

Flame




Shape No.	862	862	863	863	863	863	879	878K	879K	878	834	878
Size (1/10 mm)	012	016	012	014	016	018	018	012	012	012	021	014
Length (mm)	8.0	8.0	10.0	8.0	10.0	9.0	9.0	8.0	10.0	8.0	6.0	8.0
Extra Fine	0099	0103	0107									
Fine	0100	0104	0108	0110	0112	0114	0116			0123		
Medium	0101	0105	0109					0119	0121	0124	0173	0125
Coarse	0102	0106	0019	0020	0113	0115	0117	0120	0122	0021		0126
Super Coarse				0111			0118			0022		

Beveled



Shape No.	878
Size (1/10 mm)	012
Length (mm)	8.0
Extra Fine	
Fine	0123
Medium	0124
Coarse	0021
Super Coarse	0022

Depth Cutter



Shape No.	834
Size (1/10 mm)	021
Length (mm)	6.0
Extra Fine	
Fine	
Medium	0173
Coarse	
Super Coarse	

Cylinder



Shape No.	878	879
Size (1/10 mm)	014	012
Length (mm)	8.0	10.0
Extra Fine		
Fine		
Medium	0125	0127
Coarse	0126	0128
Super Coarse		

Bulk Reducer



Shape No.	6055
Size (1/10 mm)	018
Length (mm)	8.0
Extra Fine	
Fine	
Medium	
Coarse	0129
Super Coarse	0023

Curettage



Shape No.	878K	878K	878K	879K
Size (1/10 mm)	014	016	018	018
Length (mm)	9.0	9.0	8.0	11.0
Extra Fine				
Fine				
Medium	0130	0132	0134	0136
Coarse	0131	0133	0135	0137
Super Coarse				

Double Cone



Shape No.	806	813
Size (1/10 mm)	012	014
Length (mm)	2.5	2.0
Extra Fine		
Fine		
Medium	0138	0139
Coarse		
Super Coarse		

Single-Use Diamonds

Available in Packs of 10 & 25

Inverted Cone



Shape No.	805	805	805
Size (1/10 mm)	010	014	016
Length (mm)	1.0	2.0	2.0
Extra Fine			
Fine			
Medium	0140	0142	
Coarse	0141	0143	0144
Super Coarse			

Modified Shoulder



Shape No.	837KR	850KR
Size (1/10 mm)	012	018
Length (mm)	8.0	9.0
Extra Fine		
Fine		
Medium		
Coarse	0145	0146
Super Coarse		

Mosquito

Shape No. 392
Size (1/10 mm) 016
Length (mm) 5.0
Extra Fine
Fine 0147
Medium 0148
Coarse
Super Coarse

Pear



Shape No.	830	830
Size (1/10 mm)	010	014
Length (mm)	2.2	3.0
Extra Fine		
Fine		
Medium	0149	0151
Coarse	0150	0152
Super Coarse		0153

Round




			Large Head
Shape No.	801	801	801
Size (1/10 mm)	010	014	023
Length (mm)	1.0	1.4	2.3
Extra Fine		0156	
Fine			
Medium	0154	0157	0166
Coarse	0155	0158	
Super Coarse			

Round End Cylinder



Shape No.	838	880	880	881
Size (1/10 mm)	012	010	012	014
Length (mm)	4.0	7.0	7.0	8.0
Extra Fine				
Fine				
Medium				
Coarse	0159	0160	0162	0164
Super Coarse		0161	0163	

Wheel




	Large Head	
Shape No.	909	909
Size (1/10 mm)	035	043
Length (mm)	1.1	1.3
Extra Fine		
Fine		
Medium		
Coarse	0169	0171
Super Coarse	0170	0172

Turbo Diamonds


Available in Packs of 10

Occlusal Reduction/Finishing




Shape No. 368T
Size (1/10 mm) 023
Length (mm) 5.0
Item # Q2080174

Flat End Taper




Shape No.	847T	847T
Size (1/10 mm)	016	018
Length (mm)	9.0	8.0
Item #	Q2080175	Q2080176

Round End Taper




Shape No.	850T	850T	855T	855T	856T
Size (1/10 mm)	014	018	016	018	018
Length (mm)	10.0	10.0	7.0	7.0	9.0
Item #	Q2080177	Q2080178	Q2080179	Q2080180	Q2080181

Flame



Shape No. 863T
Size (1/10 mm) 014
Length (mm) 8.0
Item # Q2080182

Modified Shoulder



Shape No.	837KR	850KRT
Size (1/10 mm)	012	018
Length (mm)	8.0	9.0
Item #	Q2080183	Q2080184

Round End Cylinder



Shape No. 880T
Size (1/10 mm) 010
Length (mm) 7.0
Item # Q2080185

Quala® DENTAL PRODUCTS