Quala®
DENTAL PRODUCTS

RESTORATIVE

# DISPOSABLE BRUSH APPLICATORS

Bendable, disposable applicator brushes for applying all types of liquids. Assorted color handles help identify materials.

| FEATURES | BENEFITS |
|---|---|
| Quality construction | Rugged, bendable, one-piece disposable applicator brushes |
| Design | Wire-core construction with unique “easy-to-bend” design to assure retention of shaped angle |




| CODE | DESCRIPTION | QTY |
|---|---|---|
| 204002 | Regular (Black Filaments) | 144 |

**Distributed by:**

Quala®
DENTAL PRODUCTS
www.quala.com
#QSS012R3

© 2022 NDC, Inc.